



BLACK · DIAMOND

BY OSSEUS

POCT

PRECISION THREADS

Self-tapping screw tip facilitates effortless bone purchase while robust thread profile increases fixation strength.

SYSTEM VERSATILITY

Ø3.5mm and Ø4.0mm rods with standard, smooth shank, and reduction screw options allow for maximum surgical versatility.

CT-JUNCTION BRIDGE

A combination of novel connectors, unique fixation methods, and anatomic rods offers unparalleled options for bridging the cervico-thoracic junction.

ANATOMIC ANGULATION

Additional 20° of bi-lateral preferred angulation improves surgical flexibility during lateral mass screw insertion.



This Device is not available for sale in the United States

The Black Diamond POCT is a leading edge posterior occipito-cervico-thoracic screw system focused on precise screw placement, anatomic angulation, and surgical versatility. With unique options for bridging the CT-junction, the Black Diamond POCT is designed to integrate cohesively with the Black Diamond thoraco-lumbo-sacral Pedicle Screw System allowing for a seamless construct from the occiput to the sacrum.

TECHNICAL SPECIFICATIONS

BLACK DIAMOND POCT SCREWS




STYLES:	DIAMETERS:	LENGTHS:
Standard	Ø3.5mm	10mm 18mm
Smooth Shank	Ø4.0mm	12mm 20mm
	Ø4.5mm	14mm 28mm
		16mm 30mm

BLACK DIAMOND POCT RODS



MATERIALS:	DIAMETERS:	LENGTHS:
Titanium Alloy	Ø3.5mm	20mm - 120mm
Cobalt Chrome	Ø4.0mm	(5mm Increments)
	Ø5.5mm	
	Transition Rod	

BLACK DIAMOND POCT CROSSLINKS



STYLES:	SIZES:	
Rod-to-Rod	S	L
Head-to-Head	M	XL

This Device is not available for sale in the United States