

BONE PURCHASE

Aggressive distal cutting flutes and double helix threads allow for effortless bone purchase and screw insertion without pre-tapping.

FIXATION

Wide range of screw diameters and lengths offer an array of fixation options.

OSSECINTEGRATION

Proprietary surface treatment and open architecture promote bone growth onto, into, and throughout the device.

COMPRESSION


Optional lag caps facilitate compression across SI joint.



The Blue Topaz Sacroiliac Screw System consists of fenestrated, double helix threaded screws with aggressive distal cutting flutes facilitating effortless bone purchase and screw insertion. The implant's proprietary surface treatment and open architecture promote bone growth onto, into, and throughout the device. The screws come in a wide range of diameters and lengths offering an array of fixation options while attachable lag caps allow for compression across the SI joint.

TECHNICAL SPECIFICATIONS

BLUE TOPAZ SCREWS



DIAMETERS

Ø 9.5mm Ø10.5mm

LENGTHS 35mm

40mm 45mm 50mm 55mm

60mm

BLUE TOPAZ LAG CAPS



SIZES

Standard