



# WHITE PEARL™

BY OSSEUS™

## PRECISE PLACEMENT

Patented cannulated design allows for quick and accurate preliminary fixation at midline of plate.

## EXTREME ANGULATION

System framework offering  $\pm 20^\circ$  cephalad / caudal angulation allows for maximum flexibility on challenging alignments.

## ROBUST LOCKING

Manual, one-step locking rivets provide a visual, tactile, and audible confirmation that screws have been locked.

## LOW PROFILE

Ultra-thin plate designed for decreased anatomic disruption and increased disc space visibility.



